

## Comparative Analysis of N-1 Perfusion Strategies for Enhanced Protein Expression and Cost Efficiency in Biopharmaceutical Manufacturing

**Background:** Biopharmaceutical manufacturing relies on efficient cell culture strategies to achieve high yields of therapeutic proteins. Traditionally, conventional fed-batch culture methods have been widely used, but they come with challenges such as seed train scale-up requirements and long production durations. N-1 perfusion culture has emerged as an alternative strategy offering distinct approaches to address these challenges.

N-1 perfusion culture involves culturing cells at high densities (N-1 stage) by perfusion before transitioning to the production phase. This approach allows for the optimization of cell health and productivity before entering the main production phase. Within N-1 perfusion culture, there are two notable strategies:

- High Dilution During Production:** In this strategy, cells are diluted during the production phase to match the cell density observed in conventional fed-batch processes. This approach aims to achieve comparable production titers (g/L) while reducing the need for extensive seed train scale-up, thereby streamlining the overall manufacturing process.
- High Cell Density During Production:** Conversely, this strategy focuses on maintaining high cell density during the production phase. By keeping cells at higher concentrations, the production titer per unit volume can be increased significantly. This not only enhances productivity but also reduces the overall production duration, leading to improved cost efficiency.

**Methods:** The comparative analysis will involve conducting bioreactor experiments to evaluate the performance of each N-1 perfusion strategy. Key parameters such as cell viability, specific productivity, product titer, and production duration will be monitored and analyzed. Cost analysis will also be conducted to assess the economic implications of each strategy.

### Conventional 14 day-Fedbatch

- Cell thaw 50mL
- N-2: 80mL
- N-1: 300mL
- 3L (Fedbatch: 14day)

### High Dilution N-1 Perfusion

- Cell thaw 50mL
- N-1: 30ml (Perfusion)
- 3L (Fedbatch: 14day)

### High Seed N-1 Perfusion


- Cell thaw 50mL
- N-3: 80mL
- N-2: 300mL
- N-1: 1.5L (Perfusion)
- 3L (Fedbatch: 9-14day)

**Equipment:** Applikon 3L with my-Control controller, ATF2 (Repligen)

**Product A:** Symmetric Bispecific (172KD) ; **Product B:** Symmetric Bispecific (150KD)

**Cell Type:** CHO-K1 (ATCC)

**Culture Medium:** CD-CHO for seed train (Commercial Brand) ,HSM16 for N-1 perfusion and production scale (from Healsun)



**Objective:** This study aims to provide a comprehensive comparison of the two N-1 perfusion strategies: high dilution during production and high cell density during production. The comparison will be based on key performance metrics such as specific productivity, production titer, production duration, and overall process costs. This poster presentation will present a detailed comparative analysis of the two N-1 perfusion strategies, shedding light on their respective strengths and limitations. The findings from this study will contribute valuable insights to optimize N-1 perfusion culture strategies for improved bioprocess performance and economic viability in biopharmaceutical manufacturing.


### Results:

Product A	N Seed Density	Dilution Factor	Production Culture Duration (Day)	Titer (g/L)	Cell Viability During Harvest
Conventional Production	0.6-1.0E+06	7-8	14	6.10	93.40%
N-1 High Dilution Production	0.6-1.0E+06	46	14	6.11	91.00%
N-1 High Seed Production	3.60E+07	2	9	8.71	81.80%

Product B	N Seed Density	Dilution Factor	Production Culture Duration (Day)	Titer (g/L)	Cell Viability During Harvest
Conventional Production	0.6-1.0E+07	7-8	14	5.43	94.00%
N-1 High Dilution Production	0.6-1.0E+07	80	14	5.61	93.30%
N-1 High Seed Production	6.20E+07	2	9	9.06	85.80%

### Cost Estimation Analysis for 2000L Scale:

High dilution fedbatch	High seed fedbatch	Traditional fedbatch
Cell thaw 50mL	Cell thaw 50mL	Cell thaw 50mL
N-3 300mL	N-5 300mL	N-5 300mL
N-2 2L	N-4 2L	N-4 2L
N-1 20L(Perfusion)	N-3 10L	N-3 12L
2000L(fedbatch,titer:5-7g/L)	N-2 40L	N-2 70L
Time: 14day	N-1 250L	
Save equipment costs	N-11000L(Perfusion) & 2000L(fedbatch, titer:8-12g/L) time:9-14day	2000L(fedbatch, titer:5-7g/L) time:14day
Shorten cultivation time Enhance production		
Process Map for 2000L scale		



### 1. High Dilution N-1 Perfusion vs. Conventional 14 day-Fedbatch

Take 2000L scale as an example for product A culture medium cost calculation

Materials	Item	Traditional fedbatch Cost of a batch (2,000L)	High dilution fedbatch Cost of a batch (2,000L)	Item	Traditional fedbatch	High dilution perfusion
Seed culture medium (commercial)	Unit Price (USD/L)	28	28	Total Price per Batch (USD)	61,340	58,956
	Quantity (L)	80	2			
	Cost (USD)	2,240	56			
Basal medium (N-1 stage)	Unit Price (USD/L)	10	10	Expression Level (g/L)	6.10	6.11
	Quantity (L)	350	130			
	Cost (USD)	3,500	1,300			
Basal medium (FB stage)	Unit Price (USD/L)	10	10	Batch Yield (g)	12,220	12,200
	Quantity (L)	1,400	1,600			
	Cost (USD)	14,000	16,000			
Feeding medium 1	Unit Price (USD/L)	97	97	Cost per Gram of Protein (USD)	5.0	4.8
	Quantity (L)	400	400			
	Cost (USD)	38,800	38,800			
Feeding medium 2	Unit Price (USD/L)	70	70			
	Quantity (L)	40	40			
	Cost (USD)	2,800	2,800			

Cost analysis reveals compelling economic benefits associated with High dilution N-1 perfusion culture. While material costs may not significantly differ between high dilution N-1 perfusion and traditional FB processes, **considerations of labor and equipment operation costs unveil a distinct advantage for the high dilution approach.** The lower batch cost of the high dilution process, attributed to reduced labor and equipment operation expenses, positions it as a cost-effective solution. The adoption of N-1 high dilution perfusion presents a pathway to expedited processes capable of yielding multi-gram titers, thereby significantly reducing time constraints by saving seed trains and overall costs. This accelerated production not only facilitates swift pipeline delivery but also optimizes facility utilization.

### Conclusion:

Regarding N-1 high dilution perfusion, N-1 high seed perfusion, and their comparison with conventional fed-batch processes, along with the importance of selecting a suitable culture medium vendor:

#### 1. N-1 High Dilution Perfusion Conclusion:

- N-1 high dilution perfusion offers substantial benefits in terms of seed train saving
- The process allows for faster production, improved facility utilization.
- Cost analysis reveals that while material costs may not differ significantly from traditional fed-batch processes, considerations of labor and equipment operation costs demonstrate a cost-effective advantage for the high dilution approach.
- It is crucial to select a culture medium vendor that offers competitive pricing, as the cost of each gram of product can be impacted by the choice of culture medium.

#### 2. N-1 High Seed Perfusion Conclusion:

- N-1 high seed perfusion also presents notable advantages such as reduced seed train requirements, shorter production durations, and increased product titers close to twice that of conventional fed-batch processes.
- While there is higher consumption of basal medium during the seed train stage, overall cost savings can be achieved through reduced seed train time and production days.
- Similar to high dilution perfusion, the selection of a culture medium vendor is crucial in optimizing costs and ensuring economic viability.

#### 3. Comparison with Conventional Fed-Batch Processes:

- Both N-1 high dilution perfusion and N-1 high seed perfusion outperform conventional fed-batch processes in terms of efficiency, productivity, and cost-effectiveness.
- The choice between the two N-1 perfusion strategies depends on specific requirements such as seed train management, production duration, and desired product titers.
- The selection of a culture medium vendor plays a vital role in determining the overall cost per gram of product, highlighting the importance of vendor assessment and cost analysis in biopharmaceutical manufacturing.

#### 4. Scale-Up Considerations and Challenges with Larger ATF Systems:

- Scaling up N-1 perfusion culture, particularly with larger-scale ATF (Alternating Tangential Flow) systems, is essential for meeting production demands.
- However, it's important to note that as the scale of the ATF system increases, certain challenges may arise, particularly with large hollow fiber systems.
- Issues such as membrane fouling, uneven flow distribution, and potential blotting problems can occur in larger-scale hollow fiber systems, impacting overall process efficiency and productivity.
- Therefore, careful attention and optimization strategies are necessary to mitigate these challenges and ensure smooth operation of larger ATF systems in N-1 perfusion culture.

**References** Mahé, A., Martiné, A., Fagète, S. et al. Exploring the limits of conventional small-scale CHO fed-batch for accelerated on demand monoclonal antibody production. Bioprocess Biosyst Eng 45, 297–307 (2022). <https://doi.org/10.1007/s00449-021-02657-w>



healsunbd@hs-biopharm.com  
www.healsun-bio.com